Total Organic Carbon Analyzer

TOC-L

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CONTENTS

4	The Utility of Shimadzu TOC Analyzers
6	Features
7	Flow Line Diagrams
8	PC-Controlled Models
10	Standalone Models

11	Space-Saving and Energy-Saving Design, Other Functions and Features
12	ASI-L Autosampler and OCT-L 8-Port Sampler
13	TNM-L Total Nitrogen Unit and SSM-5000A Solid Sample Combustion Unit

14	Compatibility of Options and Special Accessories
15	Specifications
16	External Dimensions Diagrams

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TOC-L
TOTAL ORGANIC CARBON ANALYZER

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Items	High-Sensitivity Model		Standard Model	
Model	TOC-L _{CPH}	TOC-L _{CSH}	TOC-L _{CPN}	TOC-L _{CSN}
Measurement Method	680 °C combustion catalytic oxidation - non-dispersive infrared detection (NDIR) method			
Operation Method	PC-controlled	Standalone* ¹	PC-controlled	Standalone* ¹
Measured Items	TC, IC, TOC (= TC-IC), NPOC (TOC measurement via acidification and sparging) *Option: POC, TOC (= NPOC + POC), TN			
Applicable Samples	Aqueous (optional solid/gas samples)			
Measurement Range	TC : 0 to 30,000 mg/L, IC : 0 to 35,000 mg/L (Option) TN : 0 to 10,000 mg/L, POC : 0 to 500 mg/L		TC : 0 to 30,000 mg/L, IC : 0 to 3,000 mg/L (Option) TN : 0 to 10,000 mg/L, POC : 0 to 500 mg/L	
Detection Limit	TC, IC : 4 µg/L, TN : 5 µg/L		TC : 50 µg/L, IC : 4 µg/L, TN : 20 µg/L	
Reproducibility	TC, IC, NPOC : CV 1.5 % max. or ±4 µg/L max (Optional TN : CV 3.0 % max. or ±5 µg/L max)		TC, NPOC : CV 1.5 % max. or ±50 µg/L max, IC : CV 1.5 % max. or ±4 µg/L max (Optional TN : CV 3.0 % max. or ±20 µg/L max)	
Measuring Time	TC : approx. 3 min, IC : approx. 3 min (Optional TN : approx. 4 min)		TC : approx. 3 min, IC : approx. 4 min (Optional TN : approx. 4 min)	
Sample Injection	Automatic sample injection using a syringe pump and slide type injection mechanism			
Sample Injection Volume	10 to 2,000 µL variable		TC : 10 to 150 µL variable, IC : 10 to 450 µL variable	
IC Removal	Automatic addition of acid and sparging			
Sample Dilution	Dilution rate of 2× to 50× (automatic sample dilution by syringe pump), dilution accuracy: ±2 % max. (2× to 20×), ±5 % max. (21× to 50×)			
Display and Operations	Operated by PC	Operation by color LCD screen and keyboard* ²	Operated by PC	Operation by color LCD screen and keyboard* ²
External Memory (Standalone Type)	-	USB flash memory used	-	USB flash memory used
Printer (Standalone Type)	PC USB printer can be used			
Carrier Gas	High-purity air (CO, CO2, HC content: Each 1 ppm max., dew point : -50 °C max.) Supply pressure: 200±10 kPa (Additional use of optional carrier gas regulator : 300 to 600 kPa) Optional use of nitrogen gas (not possible in the TN measurement). With the standard model, optional use of pressurized gas.			
Gas Consumption	150 mL/min (230 to 250 mL/min during sparging) (Variable flow rate)		230 mL/min (A separate 100 mL/min is required for sparging with ASI-L. (variable flow rate))	
Power Supply	100 to 240 V AC, 600 VA (Permitted range: 90 to 264 V AC)			
Applicable Regulations	CE			
Ambient Temperature Range	5 to 35 °C			
Dimensions	W340 × D660 × H480 mm (Excluding protrusions. For details, see the External Dimensions Diagram.)			
Weight	Approx. 35 kg			

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Vial Types	Select from three types: 9 mL, 24 mL, 40 mL
Number of Vials	9 mL: 93, 24 mL: 93, 40 mL: 68
Vial Septum	With dedicated septum (excluding 9 mL vials)
Sample Sparging	Possible (The optional external sparging kit is required.)
Dimensions	W370 × D540 × H490 mm (excluding protrusions)
Weight	Approx. 14 kg

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Number of OCT-L Units Connected	Up to 2 OCT-L units for a single TOC-L
Vial Types	Any sample container can be used
Number of Vials	8 for a single OCT-L 16 for dual OCT-L
Sample Sparging	No sparging with OCT-L; sparging is done in the TOC-L syringe
Dimensions	W245 × D245 × H440 mm (excluding protrusions)
Weight	Approx. 3.5 kg

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
Measurement Method	Chemiluminescence
Measured Items	TN(total nitrogen)
Measurement Range	0 to 10,000 mg/L
Detection Limit	5 µg/L(CPH, CSH) 20 µg/L(CPN, CSN)
Reproducibility	CV 3 % max.
Measuring Time	Approx. 4 min
Ozone Source Gas	Air (compressed air or housing air) 500 mL/min
Dimensions	W270 × D240 × H160 mm (excluding protrusions)
Weight	Approx. 6 kg

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


TC Oxidation Method	Combustion catalytic oxidation (TC furnace temperature: 900 °C)
IC Reaction Method	Acidification (IC furnace temperature: 200 °C)
Measured Items	TC, IC, TOC
Measurement Range	TC : 0.1 to 30 mg carbon (1 to 20 µg carbon in high-sensitivity measurement) IC : 0.1 to 20 mg carbon
Maximum Sample Amount	1 g (aqueous content : < 0.5 g)
Measuring Time	Normally 5 to 6 minutes
Carrier Gas	99.9 % O ₂ at 500 mL/min High-purity O ₂ gas is required for high-sensitivity measurement.
Power Requirements	100 to 127 or 220 to 240 V AC as ordered, 700 VA
Dimensions	W450 × D656 × H290 mm
Weight	Approx. 30 kg

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